



DANAGENE CLEAN PCR KIT

REF.0501.1 50 PURIFICATIONS
REF.0501.2 250 PURIFICATIONS
REF.0501.3 1000 PURIFICATIONS

1. INTRODUCCION

DANAGENE CLEAN PCR KIT is a versatile and effective method **for a quick PCR products clean up**. This method allows the purification of DNA up to 15 µg.

The purification is based on the selective absorption of nucleic acids in spin mini-columns with silica membranes with chaotropic agents presence.

The purified DNA is eluted in a small buffer volume (5 mM Tris-HCl pH 8.5) and can be used in enzymatic manipulations, including sequencing, cloning, restriction analysis, ligations, labelling and *in vitro* transcription.

DANAGENE CLEAN PCR KIT allows a DNA recovering of .01-10 Kb with an efficiency that varies between 75-90% depending on the fragment size.

2. KIT COMPONENTS

	Ref. 0501.1 50 preps	Ref. 0501.2 250 preps	Ref. 0501.3 1000 preps	
Binding Solution	15 ml	60 ml	240 ml	Room temperature
Wash Solution	10 ml	50 ml	200 ml	Room temperature
Elution Solution	2 ml	10 ml	40 ml	Room temperature
Spin Columns	50 unid.	250 unid.	1000 unid.	Room temperature
Recollection Tubes	50 unid.	250 unid.	1000 unid.	Room temperature

Equipment and additional reagents required

- * Microtubes.
- * Microcentrifuge.
- * Water bath.
- * Ethanol 100 %.
- * Isopropanol.

3. PROTOCOL

The protocol implies the following steps:

- **DNA binding:** the sample that has to be purified goes through the Spin column and this promote the DNA binding to the silica membrane.
- **DNA purification:** washing with the wash solution that contains ethanol and removes contaminants, oligos, proteins, etc.
- **DNA recovering:** the DNA is eluted with 5 mM Tris-HCl pH 8.5 or sterile water pH 8.5.

3.1 Preparaciones preliminares

- **ADD 40 ml** 50 test, **200 ml** 250 test o **800 ml** 1000 test **OF 100% ETHANOL TO THE WASH SOLUTION**

Label the container and keep it closed to avoid ethanol evaporation.

- **ADD 10 ml** 50 test , **40 ml** 250 test o **160 ml** 1000 test **OF 100% ISOPROPANOL TO THE BINDING SOLUTION**

Label the container and keep it closed to avoid isopropopanol evaporation.

- Pre-heat the Elution Buffer at 70°C.

3.2 Protocol

1. **Add 4 volumes Binding Solution with isopropanol** to 1 volume of PCR (50 -100 μ l). For example, 50 μ l PCR mix + 200 μ l Binding buffer with isopropanol. Mix well by pipetting.
2. **Transfer the sample to a spin column** with your recollection tube.
3. **Centrifuge for 1 minute at 10.000-12.000 r.p.m.**
4. **Remove the filtrate and add 600 μ l of the Washing solution** . Centrifuge for 1 minute at 14.000 r.p.m.
5. **Remove the filtrate and add 200 μ l of washing solution.** Centrifuge for 1 minute at 14.000 r.p.m.
6. **Remove the residual ethanol** by centrifugation for 3 minutes at 14.000 r.p.m.
7. Transfer the spin column into a new receiver tube and add **at least 25 μ l of pre-warmed Elution Buffer at 70°C**. A second step of elution will increase performance.
8. **Incubate for 2 minutes and centrifuge for 1 minute at 14.000 r.p.m.**

4. PROBLEM GUIDE AND POSSIBLE ANSWER

For any doubts or additional questions about the protocol, please contact the technical service of DANAGEN-BIOTED S.L info@danagen.es