



### Protocol for the extraction of viral RNA from viral transport medium

- 1. Add 200 µl of sample** into a microtube.
- 2. Add 400 µl Viral Lysis Buffer.** Close the microtube and vortex vigorously for 20 seconds.
- 3. Incubate at room temperature for 10-15 minutes.**
- 4. Centrifuge at maximum speed for 2 minutes.**
- 5. Transfer the supernatant to** a new microtube avoiding touching the pellet that can form.
- 6. Add 350 µl Ethanol 100% .Mix well.**
- 7. Transfer half of the lysate to a Spin column with your collection tube. Centrifuge at 8,000 rpm for 30 seconds.**
- 8. Pass the other half and centrifuge at 8,000 rpm for 30 seconds.**
- 9. Add 100 µl of RNA Wash Buffer 1.** Centrifuge at maximum speed for 1 minute.
- 10. Add 700 µl of RNA Wash Buffer 2.** Centrifuge at maximum speed for 1 minute.
- 11. Centrifuge for 3 minutes at maximum speed (11.000 x g)** to completely dry the membrane. Place the column in a 1.5 ml microtube. (not supplied)
- 12. Elute with 50 µl Elution Buffer or nuclease-free water added to the center of the column membrane.** Incubate at room temperature for 2 minutes.  
**It is very important to add the nuclease-free water in the center of the membrane so that it is completely wetted.**
- 13. Centrifugate at 10.000 rpm for 60 seconds.** **Collect the 50 µl and re-deposit in the center of the membrane.** This increases the performance.
- 14. Incubate 2 minutes and centrifuge at maximum speed.**  
If problems are observed in subsequent detections, the volume of the eluate added to the PCR or RT-PCR can be changed.

For any question regarding the work protocols o problems. Please, contact DanaGen-BioTed technical service for any comment or question regarding the protocol [info@danagen.es](mailto:info@danagen.es)