



## Protocol for total RNA isolation and genomic DNA from preserved saliva samples

### Tota RNA isolation

1. Mix well the tube containing 1 ml of preserved saliva sample in the **DANASALIVA RNA** Sample Collection Kit. **It is important** to see a homogeneous solution.
2. In a microtube with **25 ul of Proteinase K** add **600 ul of preserved saliva**. **Incubate at 55°C for 45 minutes**.
3. Incubate at **90°C for 15 minutes**.
4. Place the microtubes at 4°C for 10 minutes.
5. Add **200 µl Precipitation Buffer**. **Vortex**.
6. **Centrifuge at 14.000 rpm for 5 minutes**.
7. Pour the **supernatant** in a new 1.5 ml microtube containing **600 µl of Isopropanol**. Avoid catching the possible superficial layer and / or pellet that can be formed.
8. **Centrifuge at 14.000 rpm for 3 minutes**. The RNA will be visible as a white pellet.
9. Remove the supernatant. Add **600 µl of Ethanol 70 %** and invert the tube several times to wash the RNA pellet.
10. **Centrifuge at 14.000 rpm for 1 minute**. Carefully remove the supernatant without touching the DNA/RNA pellet. It can be re-centrifuged briefly to collect the last drops of residual ethanol.

### Removal gDNA contaminant

1. Add **100 µl of Nuclease-free-water** to the DNA/RNA pellet and resuspend by micropipette.
2. **Incubate at 55°C for 15 minutes**.
3. **Add 400 µl of Binding Buffer + 50 µl of ethanol 100%**. Mix well.
4. **Transfer the sample to a gDNA removal column**. Place the column into a collection tube.
5. **Centrifuge for 1 minute at 10.000 r.p.m.**
6. Add **400 µl de ethanol 100%** to the supernatant collected in point 5. Mix well. **Save the gDNA removal column**.

7. Take one **RNA column and its collection tube** and add the mixture from point 6. Centrifuge at **8.000 -10.000 rpm** for **60 seconds**. Pass the sample twice as the volume exceeds the capacity of the column.
8. Add **700 µl Wash Buffer**. Centrifuge at maximum speed for 1 minute.
9. Centrifuge for **3 minutes** at **maximum speed (11.000 x g)** to completely dry the membrane. Place the column in a 1.5 ml (not supplied) microtube.
10. Add **30-50 µl of Nuclease-Free water**. Incubate for 2 minutes and centrifuge **at maximum speed for 1 minute**.

### **Genomic DNA isolation**

1. Add **700 µl etanol 80%** to **the gDNA removal column-**
2. Centrifuge at maximum speed for 1 minute.
3. **Centrifugar 3 minutos a máxima velocidad** para eliminar todo el etanol.
4. Centrifuge for **3 minutes** at **maximum speed (11.000 x g)** to completely dry the membrane. Place the column in a 1.5 ml (not supplied) microtube.
5. Add **50 µl of Nuclease-Free water**. Incubate for 2 minutes and centrifuge **at maximum speed for 1 minute**.