



DANAGENE SPIN FOOD-STOOL “Bacterias” KIT

Ref. 0608.1 50 isolations

Ref. 0608.2 250 isolations

1. INTRODUCTION

1.1 Product Description

This kit has been optimized for an efficient and fast **PCR-ready bacterial DNA** extraction (*Listeria*, *Salmonella*, *E.coli*, etc) **from pre-enrichment or enrichment culture** from different **food samples, raw materials or feces** using glass fiber membrane MicroSpin columns which selectively binds the DNA.

The process includes the centrifugation of 1 ml of pre-enrichment (peptone water) or enrichment medium for concentrating the cells. The cells are lysed during an incubation time with lysozyme (supplied with the kit). After a digestion with Proteinase K (supplied with the kit) and cleaning the lysis mixture by centrifugation and DNA binding to the glass fiber packed in the MicroSpin columns. The bounded DNA is purified by several washes to remove the potential PCR inhibitors and it is finally eluted.

1.2 Kit Components

| | 50 isolations | 250 isolations | T ^a Stock |
|---------------------------------|------------------|-------------------|----------------------|
| Lysozyme Reaction Buffer | 15 ml | 75 ml | 15-25°C |
| Lysis/ Binding Buffer | 15 ml | 75 ml | 15-25°C |
| Lysozyme * | 12 mg | 60 mg | -20°C |
| Proteinase K* | 22 mg | 105 mg | -20°C |
| Desinhibition Buffer (*) | 16.5 ml | 82.5 ml | 15-25°C |
| Wash Buffer (*) | 10 ml | 50 ml | 15-25°C |
| Elution Buffer | 10 ml | 50 ml | 15-25°C |
| MicroSpin Columns | 50 unid. | 250 unid. | 15-25°C |
| Collection tubes | 100 unid. | 500 unid. | 15-25°C |

(*) These solutions must be prepared as indicated in the Preliminary Preparations section of the protocol.

ATENCION: Both the Lysis/ Binding Buffer and the Desinhibition Buffer contain Guanidine hydrochloride which is an irritant agent, for this reason we recommend to use gloves and glasses for its manipulation.

1.3 Equipment and additional reagents required

1. Isopropanol.
2. Ethanol 100 %
3. Microcentrifuge.
4. 1.5 ml microtubes

2. PROTOCOL

2.1 Preliminary Preparations

- Dissolve the Lysozyme in **1.10 ml** (50 extractions) or in **5.50 ml** (250 extractions) of nuclease-free water and store at -20°C . It is recommended to do several aliquots to avoid many thaw/freeze cycles. At this temperature it is stable for 1 year
- Dissolve the proteinase K in **1.10 ml** (50 extractions) or in **5.20 ml** (250 extractions) of nuclease-free water and store at -20°C . It is recommended to do several aliquots to avoid many thaw/freeze cycles. At this temperature it is stable for 1 year.
- Verify that the Lysis/ Binding Buffer do not have precipitates due to the low temperatures. If necessary, dissolve heating at 37°C .
- Add **10 ml** (50 extraction) or **50 ml** (250 extractions) of Ethanol 100 % to the Desinhibition Buffer. Keep the container closed to avoid the ethanol evaporation
- Add **40 ml** (50 extractions) or **200 ml** (250 extractions) of Ethanol 100 % to the Wash Buffer. Keep the container closed to avoid the ethanol evaporation.
- Pre-heat the Elution Buffer at 70°C . For some downstream applications, concentrated DNA may be required. Elution with volumes of less than $200\ \mu\text{l}$ increases the final DNA concentration in the eluate significantly but slightly reduces overall DNA yield. For samples containing $<3\ \mu\text{g}$ of DNA, eluting the DNA in $100\ \mu\text{l}$ is recommended. For samples containing less than $1\ \mu\text{g}$ of DNA, only one elution in $50\ \mu\text{l}$ is recommended.

2.2 Protocol for DNA isolation from 1 ml of pre-enrichment or enrichment medium

1. Transfer 1-1.5 ml of pre-enrichment or enrichment medium to a 1.5 ml microtube and centrifuge **for 5 minutes at 12.000-14.000 rpm.**
2. Add **300 μl of Lysozyme Reaction Buffer + 20 μl Lysozyme.** Mix well.. Incubate at 37°C for 30 minutes.
3. Add **300 μl of Lysis/ Binding Buffer + 20 μl Proteinase K.** Mix well. Incubate at 70°C for 10 minutes.
4. Add **150 μl of Isopropanol.** Mix well and centrifuge 60 seconds at 14.000 rpm .
5. Pipette the lysate into reservoir of a combined MicroSpin Column –collection tube assembly. **Centrifuge at 8.000 rpm for 60 seconds.** Remove the collection tube.
6. Place the MicroSpin column in a new collection tube and add **500 μl of Desinhibition Buffer** to the reservoir. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
7. Add **500 μl of Wash Buffer** into reservoir of MicroSpin column. **Centrifuge at 12.000 rpm for 60 seconds.** Remove the liquid.
8. 2° Wash. Add **500 μl of Wash Buffer** into reservoir of MicroSpin column. Centrifuge at 14.000 rpm for 60 seconds. Remove the liquid.
9. **Centrifuge at maximum speed for 2 minutes to remove the residual ethanol.**
10. Remove the collection tube and insert the MicroSpin column in a 1.5 ml microtube. **Add 50-200 μl of Elution Buffer** (preheated at 70°C) into reservoir of MicroSpin column. Incubate 2 minutes.
11. **Centrifuge at maximum speed for 60 seconds.** The microtube contains now genomic DNA.

3. GUPROBLEM GUIDE AND POSSIBLE ANSWER

For any question regarding the work protocols or problems. Please, contact DanaGen-BioTed technical service for any comment or question regarding the protocol.